Gebrauchsinformation – ActiMaris® OROPHARYNX Spray

1. Produktprofil

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Instructions for use – ActiMaris® OROPHARYNX Spray

1. Product profile

ActiMaris® OROPHARYNX Spray is a hypertonic alkaline ionised seawater solution with sea salt (3.0 %) and sodium hypochlorite NaOCl (0.2 %) and is effective on a bio-physical basis. The sea salt and sodium hypochlorite have an anti-inflammatory and antimicrobial effect.

ActiMaris® OROPHARYNX Spray moisturises and cleanses in interaction with the mucous membrane. ActiMaris® OROPHARYNX Spray also supports the physiological milieu on the surface of the mucous membrane through its composition of minerals, salts and trace elements. Due to its hypertonicity/hyperosmolarity, the solution has a decongestant effect, dissolves plaque and oral/pharyngeal secretions and reduces in this way the inflammatory reaction of the mucous membrane. Thanks to its cleansing effect and dissolving of secretions, swallowing difficulties and the feeling of illness are alleviated in the mentioned indications. ActiMaris® OROPHARYNX Spray is both cell-friendly and antimicrobially efficient and fast. The sea salt and sodium hypochlorite are natural preservatives and guarantee a long stability of the ActiMaris® OROPHARYNX Spray solution.

ActiMaris® OROPHARYNX Spray contains:

- no cortisone
- no xylometazoline
- no alcohol
- no antibiotic

2. Fields of application/indications

ActiMaris® OROPHARYNX Spray can be used as a supplement and support for the prophylaxis and therapy of inflammations and infections in the area of the oral cavity, pharynx and larynx, such as:

- Inflammation of the gums (stomatitis), aphthae
- After dental surgeries
- Glossitis/mucositis of the oral cavity
- Tonsillitis/tonsillopharyngitis
- Dry mouth and throat (pharyngitis granularis sicca)
- Laryngitis posterior (with reflux oesophagitis)
- Acute and chronic laryngitis (laryngitis acuta et chronica)
- Bad breath (halitosis)
- Post-operative (after surgical interventions in the mouth/throat area)
- Prophylaxis of multi-resistant germs MRSA/VRE

3. Application instructions and dosage

Press on the pump with the mouth wide open. In this way a broad spray pattern achieves an adequate volume and wets the mucous membrane of the oral cavity wall, tongue and pharyngeal and laryngeal area. Residual liquid may be swallowed.

Adults and children from 2 years:

Recommended dosage of ActiMaris® OROPHARYNX Spray: 1 – 2 sprays, 4 – 5 times daily (approx. every 4 hours) Can be used during pregnancy and lactation.

Toddlers (under 2 years) and infants:

Depending on demand, 1 – 2 sprays, 2 – 3 times daily.

Always consult your paediatrician before using the product for infants or toddlers.

At the recommended dosage, the content of 50 ml (140 µl per spray burst) is sufficient for therapy for approx. 3 weeks (approx. 360 spray bursts). When using ActiMaris® OROPHARYNX Spray, an overdose is not possible, nor is there any systemic effect. The oropharyngeal solution disintegrates into water, salt and oxygen after contact and reaction with the mucous membrane. Also, a dependency/habituation to the oropharyngeal solution is not possible.

As an alternative respectively in addition, in order to reach deeper-seated complaints respectively the lower airways, the spray solution can be inhaled with the aid of a so-called ultrasonic nebuliser (3-5 µm aerosol size) (approx. 2-3 times a day as required). To do this, unscrew the cap carefully and fill the solution onto a suitable ultrasonic nebuliser. Please consult your doctor in advance.

4. Tissue tolerance and biocompatibility

ActiMaris® OROPHARYNX Spray is suitable for long-term (>24h) and repeated use in irritated mucous membrane processes, as tissue tolerance and biocompatibility have been proven also with repeated or multiple applications. The application is almost painless and is generally well tolerated by allergy sufferers.

5. Side effects and contraindications

So far, no side effects are known.

Contraindication in case of allergy to ingredients (not observed so far).

6. Interactions

When used correctly, no interactions are to be expected - even with necessary systemic antibiotic therapies. A local combination therapy with another oropharyngeal preparation is not necessary or not recommended due to its optimal effect for the specified indications. Please consult your doctor for safety reasons, if you want to systematically use other drugs in combination.

7. General safety instructions

Keep out of the reach of children! Please avoid direct contact with coloured fabrics, as the contained sodium hypochlorite can have an oxidising or bleaching effect! Rinse immediately with water after contact with coloured fabrics.

8. Ingredients and technical data

Water, sea salt, sodium hypochlorite.

Clear, colourless, alkaline and hypertonic sea salt solution with characteristic odour; contains 3.0 % sea salt (Sal Maris) and 0.2 % NaOCl.

9. Storage instructions, shelf life, period of use

ActiMaris® OROPHARYNX Spray should be stored at room temperature between +5° and +25° C. Protect the spray from direct sunlight and/or heat.

The shelf life is 30 months. The current expiration date is printed on the packaging.

After opening the spray, the period of use is 3 months.

ActiMaris AG Manufacturer:

ActiMaris AG
Sandgrube 29
9050 Appenzell (SWITZERLAND)
info@actimaris.com
Phone +41 71 505 75 25

- Batch designation / batch number
- Order number
- Attention, observe instructions for use
- Can be used until
- Storage instructions (temperature)
- Protect from direct sunlight
- Manufacturer

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ActiMaris AG Fabricant :