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Instructions for use – ActiMaris® wound gel

1. Product profile

ActiMaris® wound gel is made of alkaline ionised seawater with sea salt (3.0 %), sodium hypochlorite NaOCl (0.2 %), and lithium magnesium sodium silicate. Sodium hypochlorite (NaOCl) and sea salt have an eliminating effect on typical sore pathogens without damaging the tissue. The decontaminating effect against MRSA/VRE tested positive. The oxidative effect of sodium hypochlorite and hypertonic saline solution are also responsible for the neutralisation of wound odours. The alkaline pH-value of ActiMaris® wound gel supports due to the provision of OH-ions (as part of a naturally balanced redox system) additionally the release of oxygen. Sea salt contains sodium chloride as its main ingredient (NaCl = 98 %) as well as 2 % other sea salts. The increased salinity of ActiMaris® wound gel in turn ensures an optimum hyperosmotic effect so that peri-wound swelling subsides quickly.

2. Fields of application

ActiMaris® wound gel is suitable: for cleaning, humidification, decontamination, biofilm resolution, reduction of swelling and physiological debridement, to create an optimal "micro environment" for:

- acute mechanical wounds (skin lacerations, bites, cuts, abrasions, fissures and contusions) and postoperative wounds
- chronic wounds (decubitus ulcers, arterial-venous ulcers, diabetic ulcers)
- nerotic, malodorous wounds and ulcerating tumours, even with cavities (hollow spaces)
- thermal and chemical wounds (burns 1st to 3rd degree)
- entry portals of urological catheters and PEG tubes and drains
- venous ulcers

for the local treatment of the skin and mucous membrane in inflammatory and infectious processes. for the prevention of infection at the mucous membrane and skin, burns and other wound types. Additionally, ActiMaris® wound gel is well suited as a wound filler.

3. Application instructions

ActiMaris® wound gel should be applied at room temperature (not above 30 °C). Before application of the wound gel, the wound should be cleaned with ActiMaris® wound irrigation solution (Forte or Sensitive). The wound gel is then liberally applied on or in the wound and covered with an inert compress, gauze etc. The duration and frequency of application should be adapted to the individual wound conditions and should be carried out until all deposits and necroses can be easily removed and the wound is visually clean. Depending on the wound condition, a bandage change can be done for non-infected wounds according to indications and phases, every 2 to 3 days. Whether necroses and coatings can be mechanically removed is depending on the condition of the wound. In cases of difficult removal of bandages, it is recommended to hydrate the wound dressings with ActiMaris® wound irrigation solution (Forte or Sensitive) for at least 5 minutes until a smooth removal of the bandages is possible without trauma to the wound surface. A primary surgical treatment is recommended with thick, firmly adhering fibrinous or necrotic wound deposits. For cleaning wounds, ActiMaris® wound gel (see also Section 9: Interactions) can be combined with inert bandages or other wound dressings.

4. Tissue tolerance and biocompatibility

ActiMaris® wound gel is suitable for long-term (> 24h) as well as for repeated use during dressing changes, even with irritated skin and mucous membranes. The application is almost painless and is generally well tolerated by allergy sufferers. Inhibition of granulation and epithelialisation is not to be feared. The wound gel is dermatologically assessed as safe.

5. Risks following improper use

- Improper use may cause softened wound edges
- Cross-contamination between wounds or patients with incorrect handling of the tube/bottle

6. Side effects

So far, in several years of clinical experience, the following symptoms occasionally occurred:

- Slight temporary burning sensation in case of sensitive wounds

7. Contraindications

Since allergies can never be excluded, ActiMaris® wound irrigation solutions or ActiMaris® wound gel should not be used with known or suspected allergies to any of its ingredients.

8. Restrictions on use

Pregnancy and lactation:

There is no found evidence of genotoxicity, carcinogenicity or embryotoxic effects of the ingredients. Due to the lack of studies and long-term clinical experience with pregnant and breast-feeding women, ActiMaris® wound gel should be used in these cases only after careful medical benefit-risk assessment.

Infants and toddlers:

Due to insufficient experience on long-term exposure, ActiMaris® wound gel should only be applied after careful medical benefit-risk assessment.

9. Interactions

ActiMaris® wound gel acts energising and oxidising as a bio-physical compound due to the content of ROS substances (O₂ + OCl⁻), which is why a simultaneous combination with other local, active substances and/or interactive wound dressings is in principle unnecessary. Contact with non-precious metals should be avoided. The wound can be covered with modern wound dressings after applying ActiMaris® wound gel.

On the other hand, a good compatibility with modern wound dressings such as silver dressings has been shown in clinical trials. Likewise, there is ample clinical evidence of good tolerability when combined with inert, non-active bandages and wound dressings such as compresses, gauze, pads, sponges, gels, hydro fibres, alginates, hydrocolloids etc. The concomitant use of ActiMaris® wound irrigation solutions (Sensitive and Forte) with ActiMaris® wound gel is possible and useful.

10. General safety instructions

ActiMaris® wound gel is suitable for external use on wounds, skin or mucous membrane. Do not use for infusion or injection! Do not swallow! Use only closed and undamaged containers. Keep out of the reach of children! It is possible for coloured fabrics to fade after contact with ActiMaris® products. Rinse immediately with water after contact with coloured fabrics.

11. Composition and technical information

Water, lithium magnesium sodium silicate, sea salt, sodium hypochlorite. Clear, colourless, alkaline and hypertonic gel with characteristic smell; contains 3.0 % sea salt (Sal Maris) and 0.2 % NaOCl.

12. Storage specification, duration, expiry date

ActiMaris® wound gel has to be stored at room temperature between +5° and +25° C. Protect from direct sunlight. The shelf life is 30 months. The current expiry date is printed on the packaging. After opening the tube, the expiration period is 3 months; always close the tube tightly after use.

ActiMaris AG Manufacturer:

ActiMaris AG
Sandgrube 29
9050 Appenzell (SWITZERLAND)
info@actimaris.com
Tel.: +41 71 505 75 25

LOT Batch designation / batch number
REF Order number
Attention, observe instructions for use
Can be used until
Storage instructions (temperature)
Protect from direct sunlight
Manufacturer

CE REF
Adriatic Health Factory d.o.o.
Nikole Tesle 17
44000 Sisak
CROATIA

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ActiMaris AG Fabricant :